


---

# Plate RUNNER HD®

## by TROPHOS

---

Rapid High Resolution  
Cell Fluorescence Imaging  
made easy



Rat Motor Neurons, zoom in this true 8192x8192 image to believe it...




[www.trophos.com](http://www.trophos.com)

Parc Scientifique de Luminy – Case 931  
Luminy Biotech Entreprises  
13288 Marseille Cedex 9 – France  
Tel : +33 (0)4 91 82 82 82  
Fax : +33 (0)4 91 82 82 89  
[info@trophos.com](mailto:info@trophos.com)

June 2009, Ed. 4g

# the TROPHOS Plate RUNNER HD® at a glance



Rat Cortex Neurons  
8192x8192 image

## The TROPHOS Plate RUNNER HD® is a lightweight compact HTS tool

- Optimized for phenotypic HTS or research assays on fluorescent adherent cells and small organisms, in 96 or 384 well plate format
- Accurately and quickly measures cell fluorescence intensity and localization,
- Enables in-depth study of cell morphology and motility in various tests (cell viability, cytotoxicity, cell cycle, membrane potential, neurite outgrowth, angiogenesis...).
- Uniquely large field of view (8mm) covering the complete well of a 96 plate
- Single shot for the full well whatever the resolution - avoiding time consuming and error prone mosaic image reconstruction,
- Excellent resolution - from standard 1024x1024 up to 8192x8192 pixel matrix giving a 1µ/pixel resolution
- Impressive scanning speed - 2mn40 for 96x1 Megapixels images up to 25mn for 96x64 Megapixels High Resolution images
- Flexible and low cost operation
  - low power consumption (10W);
  - long life LED light source (100 000h) at three commonly used wavelength (UV 365, blue 475, green 530);
  - small footprint and light weight;
  - robust - reducing dramatically the MTBF and maintenance needs,
- Standard interfaces (driven by a standard Home Laptop PC) and open data formats (full standard tiff images and csv tabular data) – no locked in, expensive software
- Full and simple automation

## The TROPHOS Plate RUNNER HD® in Competition


See the difference over other plate imaging systems based on traditional microscopy design with their small field, small CCD of only 1 Megapixels leading to mosaic images, low speed, complex and fragile architecture, over sophisticated software and closed proprietary data formats. Our simple system is now patented in USA (# 7,372,626), Europe (App.# 01934077), France (# 2,808,888) and patent is pending in Canada (App.# 2,408,731) and Japan (App.# 2004-516495). The price is also attractive (please apply to [info@trophos.com](mailto:info@trophos.com) for a quotation).

## The TROPHOS Plate RUNNER HD® - Why


For TCO reduction (including people training costs), ease of use, speed, fully automated/full plate/full well high resolution images, interoperability with any imaging software and custom data processing, unlimited offline re-analysis, and statistically significant exhaustive quantitative results.

# INSIDE THE TROPHOS Plate RUNNER HD<sup>®</sup>


**Fig. 1 Front View**



**Fig. 2 Optical System Side View**



**Fig. 3 Optical System Top View**



# TROPHOS Plate RUNNER HD<sup>®</sup> SPECIFICATIONS

<b>Sample support</b>	96, 384 well plates
<b>Spectral Ranges and Rated Powers (*)</b>	<p>Excitation :</p> <ul style="list-style-type: none"> <li>• UV – 360 / 370 nm (50 mW)</li> <li>• Blue – 455 / 490 nm (100 mW)</li> <li>• Green – 505 / 555 nm (50 mW)</li> </ul> <p>Emission :</p> <ul style="list-style-type: none"> <li>• UV <ul style="list-style-type: none"> <li>○ High Pass filter : &gt; 455 nm</li> </ul> </li> <li>• Blue <ul style="list-style-type: none"> <li>○ Band Pass Filter : 495 / 540 nm</li> <li>○ High Pass Filter : &gt; 510 nm</li> </ul> </li> <li>• Green <ul style="list-style-type: none"> <li>○ Band Pass Filter : 570 / 620 nm</li> <li>○ High Pass Filter : &gt; 570 nm</li> </ul> </li> <li>• One open position</li> </ul> <p>Typically suited for dyes (non exhaustive list):</p> <p>DAPI, Hoechst, Propidium Iodide</p> <p>Calcein, Alexa Fluor 488, FITC, CFP, GFP, Rhodamine 110, Rhodamine123, fluo-4</p> <p>DsRed, YFP, Cy3</p> <p>PE, PE-Cy5, PerCP</p>
<b>Sensitivity</b>	30 000 FITC equivalent
<b>Optical Properties</b>	<p>Field : <math>\Phi</math> 8 mm</p> <p>Numerical Aperture : 0.2</p> <p>Depth of field : about 8.5 <math>\mu</math>m at maximum resolution of 1<math>\mu</math>m, about 20<math>\mu</math>m at 2<math>\mu</math>m resolution, about 40<math>\mu</math>m at resolution of 7.4<math>\mu</math>m.</p> <p>Optical resolution : around 1 micron</p> <p>CCD resolution : 7.4<math>\mu</math>m (1 MegaPixels), 3.7<math>\mu</math>m (4 Mpx), 1.85<math>\mu</math>m (16 Mpx), 0.925<math>\mu</math>m (64 Mpx)</p>
<b>Timing Performances</b>	<p>For 200 ms exposure, on a Pentium dual core / 2GB RAM / WXP SP2 PC</p> <p>At minimum resolution of 7.4<math>\mu</math>m (1 Megapixels image size) :</p> <p>Acquisition : 2mn40s for 96 images, 5mn30 for 2x96, 6mn50s for 3x96 images (Triple color)</p> <p>Particles analysis : 7s for 96 images (2mn38s with background substraction)</p> <p>At high resolution of 2<math>\mu</math>m (16 Megapixels image size) :</p> <p>Acquisition : 8mn for 96 images, 16mn for 2x96, 23mn for 3x96 images (Triple color)</p> <p>Particles analysis : 2mn10s for 96 images (39mn with background substraction)</p> <p>At maximum resolution of 1<math>\mu</math>m (64 Megapixels image size) :</p> <p>Acquisition : 25 mn for 96 images, 51 mn for 2x96, 75 mn for 3x96 images (Triple color)</p> <p>Particles analysis : 7mn30s for 96 images (2h25mn with background substraction)</p>
<b>Hardware</b>	<p>Light Source : Powerful light emitting diodes (Lumileds<sup>™</sup>)</p> <p>Image sensor : Diagnostic Instruments DXF1500 Monochrome cooled 14bits CCD camera, 2048x2048 "μ-scanning" CCD (i.e. CCD moved by internal piezo) allowing up to 8192x8182 resolution.</p> <p>Motorized XY stage : Prior H107, Proscan controller (repeatability <math>\pm</math> 1<math>\mu</math>m , resolution 0.1<math>\mu</math>m)</p> <p>Motorized filter wheel : Thorlabs FW102B – 6 positions</p>
<b>Software</b>	<p>3 dedicated tools compatible with Windows 9x/Me/2k/XP/Vista operating systems :</p> <ul style="list-style-type: none"> <li>• Alignment setting : Align</li> <li>• Acquisition : Goelan</li> <li>• Analysis : Tina</li> </ul> <p>Provided on demand, a free of charge unsupported SDK to allow control and command from any other vendor software supporting dll (such as Universal Imaging MetaMorph<sup>®</sup> : ask directly to Universal Imaging<sup>®</sup> if you are interested in their customized add-on).</p> <p>Images are in full standard gray tiff format with no proprietary extension, allowing analysis with any other vendor software.</p> <p>Data files are in ascii csv format, compatible with any spreadsheet processor or database engine and easily readable by any custom software.</p>
<b>PC requirements</b>	Windows 98/Me/2k/XP/SP2/Vista system, with these interfaces: 2xRS 232, 1xUSB, 1xIEEE 1394 (firewire)
<b>Power supply</b>	AC 230 V and 110 V ; Standard DC 32 V + TROPHOS ICAL LED switching & regulation unit
<b>Power consumption</b>	< 250W, when used with a laptop PC
<b>Dimensions</b>	45x40x30 cm (HxWxD)
<b>Weight</b>	~ 21 Kg

\* Contact us for any other configuration required.